**Stereo-seq 芯片P for Go Spatial  
(1cm\*1cm)**

**Stereo-seq 芯片T for Go Spatial  
(1cm\*1cm)**

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---

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**Stereo-seq Chip P (1cm\*2cm,  
2cm\*2cm, 2cm\*3cm)**

**Stereo-seq Chip T (1cm\*2cm,  
2cm\*2cm, 2cm\*3cm)**

**Operation Guide For Receiving,  
Handling And Storing**

## Upon receiving

---

Stereo-seq Chips are transported by cold chain in vacuum sealed aluminum bags.

Upon receiving, check if the aluminum bag is sealed and packaged in good condition then immediately store unopened Stereo-seq Chips at  $-20\text{ }^{\circ}\text{C}$  or  $4\text{ }^{\circ}\text{C}$ .

It can be stored up to the expiration date on the label at  $-20\text{ }^{\circ}\text{C}$  or  $4\text{ }^{\circ}\text{C}$  unopened.

## Once opened

---

Place unused Stereo-seq Chips back in the slide container and in the original aluminum bag. Push out as much air as possible then keep sealed with the desiccant package inside the aluminum bag.

Check again if the aluminum bag has been sealed tightly then store at  $-20^{\circ}\text{C}$  or  $4^{\circ}\text{C}$ . Opened Stereo-seq Chips can be stored at  $-20^{\circ}\text{C}$  or  $4^{\circ}\text{C}$  up to two weeks. Please use the chips as soon as they are opened.

Please check if these happened upon opening of Stereo-seq Chips package:

- The chip type did not match the description
- The number of chips did not match the product specification
- Visible indent or scratch marks on the chip

If any of these occurred to your product, please timely report the situation to your local BD representative or FAS.

## Handling guide

---

In order to ensure the stability of transportation, Stereo-seq Chips are firmly attached to the bottom of the transparent transportation box using double-sided adhesive.

a. When removing the chip, please open the top cover of the protection box from the opening and spread it flat on the table, as shown in the picture below;



b. Insert the back end of a pair of forceps under the chip as indicated in the picture below. Then, slowly rotate the forceps to give a certain force to slowly separate the chip from the double-sided adhesive. Be cautious to detach the chip from the box gently and avoid chip from flipping over. Make sure to not scratch the front-side of the chip surface;

## Handling guide

---

c. If the chip is attached close to the inner side of the protection box, insert front side of the blunt end forceps below the chip as indicated in the picture below and gently rotate the forceps to leverage the chip out of the box;

d. Use blunt-end forceps to grip onto the lower right corner of the chip at a distance about 2~3 mm in length and width and then remove the chip from the transportation box, as shown in the figure below:

# For product-related documents

---

Please visit:

<https://en.stomics.tech/resources/documents/list>

**Step 1: Select the type of resources on the right side of the page**

**Step 2: Click on the image to preview the document. Click 'download' then sign in to download the corresponding document.**

## Contact Us

---

**Web:** <https://en.stomics.tech>

**Email:** [info\\_global@stomics.tech](mailto:info_global@stomics.tech)